
<b>FR</b>	AIGUILLE DE PONCTION LOMBAIRE MARROW CELLUTION ENCART D'INFORMATION	<b>PT</b>	AGULHA DE ASPIRAÇÃO DA MEDULA ÓSEA MARROW CELLUTION BULA
<b>GB</b>	MARROW CELLUTION BONE MARROW ASPIRATION NEEDLE INFORMATION INSERT	<b>DK</b>	MARROW CELLUTION-NÅL TIL ASPIRATION AF KNOGLEMARV INFORMATIONS BILAG
<b>US</b>		<b>SE</b>	MARROW CELLUTION™ ASPIRATIONSNÅL FÖR BENMÄRG BIPACKSEDEL
<b>DE</b>	MARROW CELLUTION KNOCHENMARK-ASPIRATIONSNÄDEL PACKUNGSBEILAGE	<b>FI</b>	MARROW CELLUTION LUUYDIN NÄYTTEENOTTONEULA TIETOSELOSTE
<b>ES</b>	AGUJA DE ASPIRACIÓN DE MÉDULA ÓSEA MARROW CELLUTION HOJA INFORMATIVA	<b>GR</b>	ΒΕΛΩΝΑ ΑΝΑΡΡΟΦΗΣΗΣ ΜΥΕΛΟΥ ΟΣΤΩΝ MARROW CELLUTION ΕΓΧΕΙΡΙΔΙΟ ΠΛΗΡΟΦΟΡΙΩΝ
<b>IT</b>	AGO PER ASPIRAZIONE DEL MIDOLLO OSSEO MARROW CELLUTION FOGLIETTO ILLUSTRAZIVO	<b>PL</b>	IGŁA DO ASPIRACJI SZPIKU KOSTNEGO MARROW CELLUTION ULOTKA INFORMACYJNA
<b>NL</b>	BEENMERG CELLUTION BEENMERG ASPIRATIE NAALD BIJSLUITER	<b>TR</b>	MARROW CELLUTION KEMİK İLİĞİ ASİRASYON İĞNESİ BİLGİ PROSPEKTÜSÜ
		<b>HU</b>	MARROW CELLUTION CSONTVELŐSZÍVŐ INFORMÁCIÓS ADATLAP
		<b>CZ</b>	JEHLA MARROW CELLUTION K ODSÁVÁNÍ KOSTNÍ DŘENĚ INFORMAČNÍ LETÁK

## FRANÇAIS

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**ENGLISH / US**

Sterile: Contents sterile unless package is opened, damaged, or expired.

**A. INTENDED PURPOSE:** The Marrow Cellution Bone Marrow Aspiration Needle is intended for aspiration of bone marrow or autologous blood using a standard piston syringe. If the Ranfac Bone Marrow Biopsy Needle is included, it is intended for the purpose of harvesting bone marrow specimens.

**B. INTENDED USER:** The Marrow Cellution Bone Marrow Aspiration Needle is intended to be used by a licensed physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of bone marrow or autologous blood aspiration. The Ranfac Bone Marrow Biopsy Needle is intended to be used by a licensed physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of bone marrow biopsy.

**C. INTENDED PATIENT POPULATION:** The Marrow Cellution Bone Marrow Aspiration Needle is intended to be used on any patients that are suitable for bone marrow or autologous blood aspiration as determined by the physician. The Ranfac Bone Marrow Biopsy Needle is intended to be used on any patients that are suitable for bone marrow harvesting as determined by the physician.

**D. INDICATION(S):** The Marrow Cellution Bone Marrow Aspiration Needle is indicated for use for aspiration of bone marrow or autologous blood using a standard piston syringe. The Ranfac Bone Marrow Biopsy Needle is indicated for the surgical harvesting of bone marrow.

**E. CONTRAINDICATION(S):** The Marrow Cellution Bone Marrow Aspiration Needle should not be used for sternal procedures. The Ranfac Bone Marrow Biopsy Needle should not be used for sternal procedures.

**F. INTENDED LIFETIME:** The Marrow Cellution Bone Marrow Aspiration Needle and Ranfac Bone Marrow Biopsy Needle are sterile, single-use devices intended for

- Transient use (<60 minutes)

**G. INTENDED CLINICAL BENEFIT:** The Marrow Cellution Bone Marrow Aspiration Needle is intended for use for aspiration of bone marrow or autologous blood using a standard piston syringe in order to minimize the need for external manipulations of the needle, to help recover adequate quantity and better quality of specimen (lower risk of sample loss) and to minimize the likelihood of complications such as patient discomfort, pain, or injury. The Ranfac Bone Marrow Biopsy Needle is intended to harvest bone marrow specimens in order to minimize the need for external manipulations of the needle, to help recover adequate quantity and better quality of specimen (lower risk of sample loss) and to minimize the likelihood of complications such as patient discomfort, pain, or injury.

**WARNINGS/CAUTIONS:** For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Do not use if package is open or damaged and if the expiry date has been exceeded. Do not continue to use if the device gets damaged during procedure. Do not bend the needle excessively, it may result to damage/breakage of the device. After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

**Note:** These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

**H. DIRECTIONS FOR USE**

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the needle from the tray and inspect the tip of cannula and stylets for any damage or other imperfections that would prevent the proper operation of the needle.
3. Determine if the black spacer that is distal of the depth stop will be necessary to achieve the desired depth. If deeper penetration is required, slide the spacer off of the needle.
4. Following proper placing technique, locate the needle adjacent to the aspiration site.
5. Determine the desired depth of insertion and set the depth stop to that distance.
6. Using gentle, but firm pressure, while holding the mated handles of the Access Needle and pointed stylet together, the user advances the Access Needle with either a clockwise/counter-clockwise twisting motion or advancement can be by use of a bone mallet.
7. Entrance into the marrow cavity is generally detected by decreased resistance.
8. Once the Access Needle is secure in the bone, remove the pointed stylet and replace it with the blunt stylet.

9. Advance the Access Needle into the marrow cavity until it reaches the desired depth. Advance the depth stop, with or without spacer, so that the tip is touching the patient.
10. Remove the blunt stylet and insert the Aspirator Cannula.
11. Attach a syringe to the needle using the luer lock connection.
12. Apply negative pressure by quickly withdrawing the syringe plunger. Aspiration will begin.
13. Hold the body of the depth stop and turn the Access Needle handle counter-clockwise one full revolution (360°). (The Aspiration Cannula & Syringe will turn with the Access Needle.) This will expose the Aspiration Cannula tip to a different marrow area.
14. When aspiration is complete, remove the syringe from the needle, and then remove the needle.

**If a biopsy needle is included, follow steps 15-23 to acquire bone marrow sample.**

15. Remove needle guard from needle and inspect the cutting edges of both cannula and stylet for any damage or other imperfections that would prevent the proper operation of the needle.
16. Following proper placing technique, locate the needle adjacent to the biopsy site.
17. Using gentle, but firm, pressure advance the needle by rotating it in an alternating clockwise-counter-clockwise motion. Entrance into the marrow cavity is generally detected by decreased resistance.
18. Unlock the stylet cap by rotating it counter-clockwise and remove the stylet from the cannula.
19. Very slowly and gently, advance the needle using clockwise-counter-clockwise motion until adequate marrow is obtained.
20. Insert the marked probe without force to check sample length in needle lumen. Remove the probe.
21. Insert the extraction cannula fully into the needle cannula.
22. Slowly rotate and remove the extraction cannula and needle cannula together with the same clockwise-counter-clockwise motion.
23. Remove the extraction cannula from the needle cannula.
24. Place the probe guide onto the distal tip of the cannula. Insert the probe into the distal end of the probe guide. Advance the probe into the needle and collect the biopsy specimen as it emerges from the opening in the handle.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**Contents:**

1 ea Needle Assembly, MC-RAN	1 ea Needle Assembly, R.J.N.*
1 ea 10cc, Luer-Lock Syringe	1 ea Probe, Molded*
1 ea Stylet, Molded, Blunt	1 ea Trap Cannula*
1 ea Aspiration Cannula	1 ea Probe Guide*

\*Contents included with MC-RAN-8C only

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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MedEnvoy UK Limited