


CERVOS KEYPRP Double (40ml)

Double Spin (CLASSIC PRP)

Materials & Preparation

CERVOS KEYPRP Double (40ml)
Incl. 2x CERVOS Process Disposable (CPD)



Additional Materials Required:

- 2x 50ml Syringe (Blood Collection / Transfer)
- 1 x 30ml Syringe (Cloudy Plasma Collection)
- 1x 10ml Syringe (PRP Resuspension / Collection)
- Phlebotomy Materials:
Butterfly (19G), Tourniquet, Bandage
- Anticoagulant i.e., ACD-A Vial or CitraFlow™ 5ml Syringe

Anticoagulant (AC) Priming Instructions (i.e., CitraFlow)

- Inject AC (5ml) to CPD-1, Flush to Prime Kit/Valve
- Extract AC from CDP-1 with 50ml Blood Collection Syringe
- Inject AC to CPD-2, Flush to Prime Kit/Valve
- Extract AC back into 50ml Blood Collection Syringe (~3ml)
- Collect Blood (~37ml)

Anti-Coagulant, Blood Collection / Processing Syringes Available Separately.

Process Steps

CPD-1

Anticoagulated Whole Blood

STEP 1
Load

- Inject Anticoagulated Blood to CPD-1, Flush to Prime Kit/Valve

Volume	A.C.	Blood
40ml	4ml	36ml

- Load CPD-1 in Centrifuge with opposing Balance Weight (equal weight)

CPD-1

Processed Blood Components

STEP 2
Process / Extract

- Process: **No Brake!**

Volume	RCF	Time
40ml	2300g	120 sec

- Remove CPD-1 from centrifuge & attach 30ml Syringe / ♀ Hub
- Extract** Plasma stopping at RBC interface (i.e., flash of RBC enters Luer Valve)

CPD-2

Consolidated Plasma

STEP 3
Load / Process

- Transfer Plasma to CPD-2*
- Adjust Balance Weight volume
- Load CPD-2 in Centrifuge with opposing Balance Weight (equal weight)
- Process: **No Brake!** 2300g for 6 Minutes

CPD-2

PPP
PRP

STEP 4
Extract PPP

- Attach 30/50ml Syringe
- Extract** Plasma until desired volume remains in CPD-2:
 - Top Line: ~2ml
 - Middle Line: ~5ml
 - Bottom Line: ~10ml

CPD-2

PRP

STEP 5
Extract PRP

- Attach 10ml Syringe
- Resuspend remaining Plasma and Platelets to uniform Suspension
- Extract PRP

Note: Swipe Luer Port with Alcohol Wipe Prior to Each Access. Do not Shake or Tilt Filled CPDs to Avoid Inadequate Mixing of Fluids.

* Ensure Luer Connection is Tightly Fastened after Removing Syringe from Process Disposable.

** Use Syringe Aspiration Technique or included Screw Set