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ENGLISH / US

Sterile: Contents sterile unless package is opened, damaged or expired.

A. INTENDED PURPOSE: The Bone Marrow Aspiration Needle is intended for use in aspirating a bone marrow specimen.

B. INTENDED USER: The Cervos Bone Marrow Aspiration Needle is intended to be used by a licensed physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of bone marrow aspiration.

C. INTENDED PATIENT POPULATION: The Cervos Bone Marrow Aspiration Needle is intended to be used on any patients that are suitable for bone marrow aspiration as determined by the physician.

D. INDICATION(S): The Cervos Bone Marrow Aspiration Needle is indicated for the use in aspirating a bone marrow specimen.

E. CONTRAINDICATION(S): The Cervos Bone Marrow Aspiration Needle should not be used for sternal procedures.

F. INTENDED LIFETIME: The Cervos Bone Marrow Aspiration Needle is a sterile, single-use device intended for

- Transient use (<60 minutes)

G. INTENDED CLINICAL BENEFIT: The Cervos Bone Marrow Aspiration Needle is intended for use in aspirating a bone marrow specimen in order to minimize the need for external manipulations of the needle, to help recover adequate quantity and better quality of specimen (lower risk of sample loss) and to minimize the likelihood of complications such as patient discomfort, pain, or injury.

WARNINGS/CAUTIONS: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Do not use if package is open or damaged and if the expiry date has been exceeded. Do not continue to use if the device gets damaged during procedure. Do not bend the needle excessively, it may result to damage/breakage of the device. After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner

is responsible for the proper procedure and techniques to be used with this device.

H. DIRECTIONS FOR USE

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the needle from the pouch and inspect the cutting edge of both cannula and stylet for any damage or other imperfections that would prevent the proper operation of the needle.
3. Following proper placing technique, locate the needle adjacent to the aspiration site.
4. Using a surgical drill with appropriate attachments, drill clockwise to advance needle. Entrance into the marrow cavity is generally detected by decreased resistance.
5. Rotate the stylet to the center position and pull the stylet out of the cannula.
6. If a blunt stylet is included, replace sharp stylet with the blunt stylet and advance the needle to the desired depth. If no blunt stylet is included, skip to Step 9.
7. Remove the blunt stylet.
8. Attach a male luer-lock syringe to the needle using the luer lock connection.
Note: Avoid using glass syringes because they may result in an inadequate seal and may break or fragment.
9. Apply negative pressure by quickly withdrawing the syringe plunger. Disengage the syringe and remove the aspirated sample.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Ranfac® Corp.
MedEnvoy UK Limited