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ENGLISH / US

Sterile: Contents sterile unless package is opened, damaged, or expired.

A. INDICATIONS FOR USE

The Bone Marrow Aspiration Needle is intended for use for aspiration of bone marrow or autologous blood using a standard piston syringe.

B. CONTRAINDICATIONS

Use only for bone marrow aspiration as determined by a licensed physician. The device is intended to be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone marrow aspiration. The procedure should be performed on patients that are suitable for such procedure only.

CAUTION: For single patient use only. The design of this device may not perform as intended by the manufacturer if it is re-used. The manufacturer cannot guarantee the performance, safety and reliability of a reprocessed device.

After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose of in accordance with applicable laws and regulations. Carefully place the used needle in a sharps biohazard container after the procedure is completed.

Note: *These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.*

C. DIRECTIONS FOR USE

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the needle from the pouch and inspect the tip of cannula and stylet for any damage or other imperfections that would prevent the proper operation of the needle.
3. Following proper placing technique, locate the needle adjacent to the aspiration site.
4. Using a surgical drill with appropriate attachments, drill clockwise to advance needle. Entrance into the marrow cavity is generally detected by decreased resistance.
5. Unlock the stylet handle by rotating it to the center position and remove the stylet from the cannula.
6. If a blunt stylet is included proceed with step 7. If there is no blunt stylet proceed to step 9.
7. Install the blunt stylet and lock it into place by rotating it clockwise. Using a surgical drill with appropriate attachments, drill clockwise to advance the needle into the marrow cavity until it reaches the desired depth.
8. Unlock the blunt stylet handle by rotating it to the center position and remove the stylet from the cannula.
9. Install the slip fit handle by aligning the large indicator arrow with the front fully slotted face of the cannula handle and sliding it down into place.
10. Advance the molded grip by unscrewing the grip in a clockwise direction while holding the needle handle until the outer cannula is secure against the bone.
11. Insert the Aspirator Cannula and lock into place using the luer lock connection. Ensure the white indicator on the aspiration cannula lines up with the fully slotted face of the cannula handle and the large indicator arrow on the slip fit handle. This will align the needle cannula holes with the aspiration cannula holes.
12. Attach a syringe to the aspiration cannula using the luer lock connection.
13. Apply negative pressure by quickly withdrawing the syringe plunger. Aspiration will begin.
14. Hold the needle handle and turn the molded grip clockwise one full revolution (360°). This will expose the Aspiration Cannula tip to a different marrow area.
15. When aspiration is complete, remove the syringe from the Aspiration Cannula, remove the Aspiration Cannula, remove the slip fit handle, install the drillable stylet and then drill in reverse to remove the needle assembly.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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Ranfac® Corp.
MedEnvoy UK Limited