

ELMI CM-7S Plus

Swing-Out Centrifuge



FEATURES



Cold room and refrigerator safe



Compact, Lightweight, and Versatile



Rotor imbalance technology



Automatic lid lock for convenience and safety



Rapid acceleration/braking system



Multiple brake speeds, 30 seconds to 2.5 minutes



Digital timer



LCD Control Panel

Control speed in RPM or G-Force

SPECIFICATIONS

Rotor Rotation Speed	100 - 3,500 RPM
Timer	1 second - 12 hours
Maximum Centrifugal Force	2,300 RCF (G-Force)
Brake Levels	6
Working Temperature Range	from +10 to +40 °C
Relative Air Moisture at 20 °C (68 °F)	80%
Voltage and Frequency	100 - 240 V; 50/60 Hz
Maximal Consumable Power	250 W
Size (L x W x H)	431 x 410 x 247 mm 17 x 16x 10"
Weight	13.2 kg/29 lbs

ACCESSORIES

SKU	Product Name	MSRP
CM-6MT-8	CM-7S PLUS Benchtop Swing-Out Centrifuge, 3500 RPM. Stainless Steel	
CM-6MT-3S	Rotor 6M.05 - 12 x 15ml Test Tubes - Steel adapters included	
CM-6MT-5S	Rotor 6M.02 - 24 x 15ml Test Tubes - Steel adapters included	
CM-6MT-7	Spacer for microtainers (3.5-10ml)-Compatible with 6M.05 & 6M.02	
CM-6MT-20	Rubber spacer for glass tubes-Compatible with 6M.05 & 6M.02	
CM-6MT-4	Rotor 6M.01 - 4 x 50ml Test Tubes - Plastic adapters included	
CM-6MT-12	Rotor 6M.06 - 6 x 50ml Test Tubes - Plastic adapters included	
CM-6MT-15	Adapter/Insert for 15ml tubes - for use with 50ml buckets from rotors 6M.01&6M.06	
CM-6MT-21	Combo Rotor 6M.07 2x50ml + 6x15ml - Adapters Included	
CM-6MT-22	Rotor 6M.08 4 x 60ml Test Tubes - Stainless steel adapters included	

ELMI North America

Toll Free: **(800)778-3570**

Email: info@elminorthamerica.com

elminorthamerica.com