

Hettich - Rotofix 32A
Rotor 1624, Cup 1481 (r = 149mm)

Radius in mm: 149

RPM	RCF (g)
500	42
1000	167
1500	375
2000	666
2500	1.041
2600	1.126
2700	1.214
2800	1.306
2900	1.401
3000	1.499
3100	1.601
3200	1.706
3300	1.814
3400	1.926
3500	2.041
3600	2.159
3700	2.281
3800	2.405
3900	2.534
4000	2.665

RCF (g)	RPM
100	775
500	1.732
1000	2.450
1500	3.001
1600	3.099
1700	3.195
1800	3.287
1900	3.377
2000	3.465
2100	3.551
2200	3.634
2300	3.716
2400	3.796
2500	3.874
2600	3.951
2665	4.000

Drucker - Boost 4+ Flex
Boost Bucket 03-1-0007-0123K

Radius in mm: 164

RPM	RCF (g)
500	46
1000	183
1500	413
2000	733
2500	1.146
2600	1.239
2700	1.337
2800	1.437
2900	1.542
3000	1.650
3100	1.762
3200	1.878
3300	1.997
3400	2.120
3500	2.246
3600	2.376
3700	2.510
3800	2.648
3900	2.789
4000	2.934

RCF (g)	RPM
100	739
500	1.651
1000	2.335
1500	2.860
1600	2.954
1700	3.045
1800	3.133
1900	3.219
2000	3.303
2100	3.384
2200	3.464
2300	3.542
2400	3.618
2500	3.693
2600	3.766
2700	3.837
2800	3.908
2900	3.977
2934	4.000

Drucker - Boost 2+ Flex
Boost Bucket 03-1-0007-0123K

Radius in mm: 139

RPM	RCF (g)
500	39
1000	155
1500	350
2000	622
2500	971
2600	1.051
2700	1.133
2800	1.218
2900	1.307
3000	1.399
3100	1.493
3200	1.591
3300	1.692
3400	1.796
3500	1.904
3600	2.014
3700	2.127
3800	2.244
3900	2.364
4000	2.486
4400	3.009

RCF (g)	RPM
100	802
500	1.794
1000	2.537
1500	3.107
1600	3.209
1700	3.307
1800	3.403
1900	3.497
2000	3.587
2100	3.676
2200	3.763
2300	3.847
2400	3.930
2500	4.011
2600	4.090
2700	4.168
2800	4.245
2900	4.320
3000	4.394
3009	1.407

ELMI CM-8S
Rotor 8M.03

Radius in mm: 165

RPM	RCF (g)
500	46
1000	184
1500	415
2000	738
2500	1.153
2600	1.247
2700	1.345
2800	1.446
2900	1.551
3000	1.660
3100	1.773
3200	1.889
3300	2.009
3400	2.132
3500	2.260
3600	2.391
3700	2.525
3800	2.664
3900	2.806
4000	2.952
4500	3.736

RCF (g)	RPM
100	736
500	1.646
1000	2.328
1500	2.852
1600	2.945
1700	3.036
1800	3.124
1900	3.209
2000	3.293
2100	3.374
2200	3.453
2300	3.531
2400	3.607
2500	3.681
2600	3.754
2700	3.826
2800	3.896
2900	3.965
3000	4.033
3500	4.356
3736	4.500

Eppendorf 5702
Rotor A-4-38 w round buckets

Radius in mm: 139

RPM	RCF (g)
500	39
1000	155
1500	350
2000	622
2500	971
2600	1.051
2700	1.133
2800	1.218
2900	1.307
3000	1.399
3100	1.493
3200	1.591
3300	1.692
3400	1.796
3500	1.904
3600	2.014
3700	2.127
3800	2.244
3900	2.364
4000	2.486
4400	3.009

RCF (g)	RPM
100	802
500	1.794
1000	2.537
1500	3.107
1600	3.209
1700	3.307
1800	3.403
1900	3.497
2000	3.587
2100	3.676
2200	3.763
2300	3.847
2400	3.930
2500	4.011
2600	4.090
2700	4.168
2800	4.245
2900	4.320
3000	4.394
3009	1.407

Hanil - Fleta 4
Rotor S85-4, buckets B85

Radius in mm: 172,6

RPM	RCF (g)
500	48
1000	193
1500	434
2000	772
2500	1.206
2600	1.304
2700	1.407
2800	1.513
2900	1.623
3000	1.737
3100	1.854
3200	1.976
3300	2.101
3400	2.231
3500	2.364
3600	2.501
3700	2.642
3800	2.786
3900	2.935
4000	3.087

RCF (g)	RPM
100	720
500	1.610
1000	2.276
1500	2.788
1600	2.880
1700	2.968
1800	3.054
1900	3.138
2000	3.219
2100	3.299
2200	3.377
2300	3.452
2400	3.527
2500	3.599
2600	3.671
2700	3.741
2800	3.809
2900	3.877
3000	3.943
3009	1.747

